

A prototype electrolysis cell has manufactured as per the customer's detail. The prototype cell was used for generation of chlorine gas for cobalt oxidation.



Project: sxk-127

Client: Metric PC Systems

Location: Texas, USA

Capacity: 100 mL per minute

Cathodes: one

Anodes: one

Rectifier: not provided

Items Provided:

- Prototype EW cell with anode and cathode compartment
- Membrane to keep anolyte separated from catholyte

SX Kinetics, Inc.

39 A Hamilton Ave.
Cobourg, Ontario, Canada, K9A 1V7

Web Site: <http://www.sxkinetics.com>

Further Information Contact:

Mr. Steven Webster
Hydrometallurgist

Tel +1 905 372 5981

Email: steven.webster@sxkinetics.com